## Park Systems Atomic Force Microscopy



# IMAGE GALLERY

Here, at Park Systems, we offer a full range of advanced imaging solutions for a wide variety of research applications. Enjoy the images in the gallery which highlight examples from a wide variety of sample types and imaging modes.



## 02

# 01. Topography

# 02. Advanced mode

C <sub>36</sub> H <sub>74</sub> on HOPG	04/05
Molecular network on HOPG	06
Isotactic polypropylene	07
Graphene/hBN heterostructure	08/09/10
tBN+MLG	11
Di-Phe-Phe nanotubes	12/13
Block copolymer	14
SBS	15
PS-PVAc	16
GeTe	17

NiO on ITO glass	— 18
Organosilane SAM	— 19
BaTiO₃ thin film	20
Defects on Si wafer	21
GaN LED wire	22
Line/Space patterns	23
Auto stitched WLI image	24/ 25
SOI wafer	26
Epitaxial gallium nitride (epi-GaN) film –	27
Dendrimer	28
Blood cell	29

MLG-hBN , tBG+bBN	30
LiNbO <sub>3</sub>	31
All-solid-state Li ion battery	32/33
PZT thin film	34
Aluminium TX630 alloy	35
CAM on HOPG	36

F <sub>14</sub> H <sub>20</sub>	- 37/38
Boron nitride on monolayer graphene –	39

SRAM	40/41
Ta/ NiFe /Ta microman	42
PS-PMMA block copolymer	43

# C<sub>36</sub>H<sub>74</sub> on HOPG

## C<sub>36</sub>H<sub>74</sub> on HOPG (250 nm scan)



#### True Non-contact<sup>™</sup> Mode

In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.

#### Height (2 µm × 1 µm scan)



#### Line profile



#### Height of position 2



#### Zoom in (250 nm × 250 nm scan)





Peak to valley: 3.0 nm RMS roughness: 0.37 nm

#### Scanning conditions

System: FX40 Scan Size: 250 nm × 250 nm Scan Mode: Non-contact Scan Rate: 0.3 Hz

## Multi-line profiles



## Molecular network on HOPG (400 & 100 & 25 nm scan)



#### Tapping Mode

In this alternative technique to non-contact mode, the cantilever again oscillates just above the surface, but at a much higher amplitude of oscillation. The bigger oscillation makes the deflection signal large enough for the control circuit, and hence an easier control for topography feedback. It produces modest AFM results but blunts the tip's sharpness at a higher rate, ultimately speeding up the loss of its imaging resolution.

## Isotactic polypropylene (500 & 100 nm scan)



## **Tapping Mode**

#### Height (500 nm scan)



Moiré pattern



Molecular lattice structure

When co-adsorbed on an atomically flat surface, melamine and cyanuric acid form 2D-molecular network with a period of 0.98 nm. Due to a lattice mismatch between molecular network and underlying HOPG lattice, a Moiré pattern can be observed as well.







Oriented isotactic polypropylene film displaying lamellar structure. Individual molecules can be seen as well.

#### Scanning conditions

System: NX20 Scan Size: Left 400 nm × 400 nm Middle 100 nm × 100 nm Right 25 nm × 25 nm

Scan Mode: Tapping Scan Rate: Left 2 Hz Middle 3 Hz Right 6 Hz

Cantilever: Multi75Al-G (k=3 N/m, f=75 kHz) Pixel Size: All 512×512

#### **Scanning conditions**

System: NX20 Scan Size: Left 500 nm × 500 nm Middle/Right 100 nm × 100 nm Scan Mode: Tapping Scan Rate: Left 2 Hz

In this alternative technique to non-contact mode, the cantilever again oscillates just above the surface, but at a much higher amplitude of oscillation. The bigger oscillation makes the deflection signal large enough for the control circuit, and hence an easier control for topography feedback. It produces modest AFM results but blunts the tip's sharpness at a higher rate, ultimately speeding up the loss of its imaging resolution.

## Height (100 nm scan)



Peak to valley: 4.1 nm RMS roughness: 0.55 nm

#### Phase (100 nm scan)





## **Graphene/hBN heterostructure**

## Graphene/hBN heterostructure (350 & 250 nm scan)



#### Tapping Mode

In this alternative technique to non-contact mode, the cantilever again oscillates just above the surface, but at a much higher amplitude of oscillation. The bigger oscillation makes the deflection signal large enough for the control circuit, and hence an easier control for topography feedback. It produces modest AFM results but blunts the tip's sharpness at a higher rate, ultimately speeding up the loss of its imaging resolution.

Large scale height image showing typical landscape of microstamped graphene layer with bubbles and wrinkles. A set of high-resolution images shows various Moiré patterns originating from lattice mismatch.





Imaging Moiré pattern on graphene in Tapping modes.

#### **Scanning conditions**

System: NX20 Scan Size: Left 350 nm × 350 nm Middle 350 nm × 350 nm Right 250 nm × 250 nm

Scan Mode: Tapping Scan Rate: Left 2 Hz Middle 2 Hz Right 3 Hz

#### Phase (250 nm scan)

## Graphene/hBN heterostructure (500 & 250 & 100 nm scan)



#### **Contact Mode**

In this method, the cantilever scans across a sample surface. Because the cantilever is in contact with the surface, strong repulsive force causes the cantilever to deflect as it passes over topographical features.

# Boron nitride on monolayer graphene (80 nm scan)



In this alternative technique to non-contact mode, the cantilever again oscillates just above the surface, but at a much higher amplitude of oscillation. The bigger oscillation makes the deflection signal large enough for the control circuit, and hence an easier control for topography feedback. It produces modest AFM results but blunts the tip's sharpness at a higher rate, ultimately speeding up the loss of its imaging resolution.





4.7 nm striated superlattice on boron nitride on top of monolayer graphene (tBN+MLG) was visualized in both height and phase signals.

• Sample courtesy: Qiong Ma, Boston College, US

#### **Scanning conditions**

System: NX10	
--------------	--

## Scan Mode: Tapping Scan Rate: 4 Hz

# Lateral (250 nm scan) Lateral (500 nm scan) Lateral (100 nm scan) 100 nm 50 nm

Imaging Moiré pattern on graphene in LFM mode based on Contact mode. Scan angle of LFM image is set to 45°.

#### **Scanning conditions**

System: NX20 Scan Size: Left 500 nm × 500 nm Middle 250 nm × 250 nm Right 100 nm × 100 nm

Scan Mode: Contact Scan Rate: Left 2 Hz Middle 3 Hz Right 3 Hz

Cantilever: Multi75Al-G (k=3 N/m, f=75 kHz) Pixel Size: All  $1024 \times 1024$ 

Scan Size: 80 nm × 80 nm

#### **Tapping Mode**

## Di-Phe-Phe nanotubes (70 & 30 & 25 nm scan)



Phase (70 nm scan)

20 nm

#### **Tapping Mode**

Phase (30 nm scan)

In this alternative technique to non-contact mode, the cantilever again oscillates just above the surface, but at a much higher amplitude of oscillation. The bigger oscillation makes the deflection signal large enough for the control circuit, and hence an easier control for topography feedback. It produces modest AFM results but blunts the tip's sharpness at a higher rate, ultimately speeding up the loss of its imaging resolution.

# **Di-Phe-Phe nanotubes (NTs)**



#### L-Phenylalanyl-L-phenylalanine, Di-L-phenylalanine





5 nm



5 nm

Phase (25 nm scan)

Height



#### **Scanning conditions**

System: NX20 Scan Size: Left 70 nm × 70 nm Middle 30 nm × 30 nm Right 25 nm × 25 nm

Scan Mode: Tapping Scan Rate: Left 3 Hz Middle 3 Hz Right 2 Hz

Cantilever: Multi75Al-G (k=3 N/m, f=75 kHz) Pixel Size: All 512 × 512

#### **Tapping Mode**

In this alternative technique to non-contact mode, the cantilever again oscillates just above the surface, but at a much higher amplitude of oscillation. The bigger oscillation makes the deflection signal large enough for the control circuit, and hence an easier control for topography feedback. It produces modest AFM results but blunts the tip's sharpness at a higher rate, ultimately speeding up the loss of its imaging resolution.



#### Large scale image of peptide NTs on mica

## Young's modulus



## **Block copolymer**



#### Tapping Mode

In this alternative technique to non-contact mode, the cantilever again oscillates just above the surface, but at a much higher amplitude of oscillation. The bigger oscillation makes the deflection signal large enough for the control circuit, and hence an easier control for topography feedback. It produces modest AFM results but blunts the tip's sharpness at a higher rate, ultimately speeding up the loss of its imaging resolution.

# **Poly(styrene-block-butadiene styrene) (SBS)**







Peak to valley: 3.97 nm RMS roughness: 0.44 nm

#### **Scanning conditions**

System: NX10 Scan Size: 2 µm × 2 µm Scan Mode: Tapping Scan Rate: 1 Hz

Cantilever: AC160TS (k=26 N/m, f=300 kHz) Pixel Size: 1024 × 512

0.8

1.2

μm

#### **Scanning conditions**

System: FX40 Scan Size: 1 µm × 1 µm Scan Mode: Tapping Scan Rate: 1 Hz

#### **Tapping Mode**

In this alternative technique to non-contact mode, the cantilever again oscillates just above the surface, but at a much higher amplitude of oscillation. The bigger oscillation makes the deflection signal large enough for the control circuit, and hence an easier control for topography feedback. It produces modest AFM results but blunts the tip's sharpness at a higher rate, ultimately speeding up the loss of its imaging resolution.

# **Polystyrene-polyvinyl acetate (PS-PVAc)**



#### Tapping Mode

In this alternative technique to non-contact mode, the cantilever again oscillates just above the surface, but at a much higher amplitude of oscillation. The bigger oscillation makes the deflection signal large enough for the control circuit, and hence an easier control for topography feedback. It produces modest AFM results but blunts the tip's sharpness at a higher rate, ultimately speeding up the loss of its imaging resolution.

Phase

25

# **Germanium telluride (GeTe)**



In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.





Peak to valley: 115.5 nm RMS roughness: 8.51 nm

Blended polystyrene and polyvinyl acetate (PS-PVAc) film.





Peak to valley: 39.9 nm RMS roughness: 3.4 nm

• Sample courtesy: Raffaella Calarco, Paul-Drude-Institut für Festkörperelektronik (PDI), Berlin

#### **Scanning conditions**

System: NX10 Scan Size: 4  $\mu$ m imes 4  $\mu$ m

Scan Mode: Non-contact Scan Rate: 0.5 Hz

#### **Scanning conditions**

System: FX40 Scan Size: 10 µm × 10 µm Scan Mode: Tapping Scan Rate: 1 Hz

Cantilever: AC160TS (k=26 N/m, f=300 kHz) Pixel Size: 512 × 256

#### True Non-contact<sup>™</sup> Mode



#### Triangular domains formed by germanium telluride when grown on silicon 7x7 by molecular beam epitaxy.

# Nickel oxide (NiO) on indium tin oxide (ITO) glass



#### True Non-contact<sup>™</sup> Mode

In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.

# **Organosilane self-assembled monolayer (SAM)**



In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.



Peak to valley: 36.0 nm RMS roughness: 7.51 nm

Surface morphology comparison of NiO coated on ITO glass by two different preparations.



#### **Scanning conditions**

System: NX10 Scan Size: 5  $\mu$ m imes 5  $\mu$ m

Scan Mode: Non-contact Scan Rate: 0.5 Hz

Cantilever: AC160TS (k=26 N/m, f=300 kHz) Pixel Size: 512 × 512

Peak to valley: 17.0 nm

RMS roughness: 2.38 nm

#### **Scanning conditions**

System: NX10 Scan Size: 1  $\mu$ m imes 1  $\mu$ m Scan Mode: Non-contact Scan Rate: 0.3 Hz

#### True Non-contact<sup>™</sup> Mode



## Barium titanate (BaTiO<sub>3</sub>, BTO) thin film



#### True Non-contact<sup>™</sup> Mode

In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.

## **Defects on Si wafer**



In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.

#### Height (3 µm scan)



Peak to valley: 32.17 nm RMS roughness: 1.27 nm

#### Height (1 µm scan)



Peak to valley: 13.56 nm RMS roughness: 1.00 nm

Defect 1 (4 µm scan)



Peak to valley: 15.34 nm RMS roughness: 1.99 nm

#### **Scanning conditions**

System: NX-Wafer Scan Size: Left 3 µm × 3 µm Right 1 µm × 1 µm Scan Mode: Non-contact Scan Rate: Left 0.5 Hz Right 1 Hz

Cantilever: AC160TS (k=26 N/m, f=300 kHz) Pixel Size: Left 2048 × 256 Right 512 × 512

#### Scanning conditions

System: NX-Wafer Scan Size: Left 4 µm × 4 µm Right 20 µm × 20 µm Scan Mode: Non-contact Scan Rate: Left 1 Hz Right 0.2 Hz

#### True Non-contact<sup>™</sup> Mode



#### Defect 2 (20 µm scan)

Peak to valley: 33.56 nm RMS roughness: 2.47 nm

## Gallium nitride (GaN) LED wire

# 10000

#### True Non-contact<sup>™</sup> Mode

In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.

## Line/Space patterns

Narrow Trench Mode



Peak to valley: 517.4 nm RMS roughness: 184.7 nm

3D

X : Y : Z scale = 1 : 1 : 1





#### **Scanning conditions**

System: NX-HDM Scan Size: 2 µm × 2 µm Scan Mode: Non-contact Scan Rate: 1 Hz

Cantilever: OMCL-AC55TS (k=85 N/m, f=1.6 MHz) Pixel Size: 2048 × 512

#### Scanning conditions

System: NX-Wafer Scan Size: 400 nm × 200 nm Scan Mode: NTM Scan Rate: 1 Hz

23

Line/Space patterns measured in Narrow trench mode (NTM).

# Auto stitched WLI image



#### WLI image

- Auto stitched
- Lens magnification: ×2.5
- Field of view: 22,510 μm × 23,930 μm





X : Y : Z scale = 1 : 1 : 8000

AFM and WLI measurements on VLSI standard, SHS-1800 QC (Chrome-coated) sample with certified step height  $183.9 \pm 2.0$  nm.

## Silicon on insulator (SOI) wafer



#### True Non-contact<sup>™</sup> Mode

In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.

# Epitaxial gallium nitride (epi-GaN) film



In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.



## Line profile



#### **Statistics**

Rpv(nm)	Rq(nm)	Ra(nm)
2.180	0.232	0.185

ov: Peak to valley : RMS roughness : Average roughness





Peak to valley: 1.2 nm RMS roughness: 0.14 nm

Extremely flat surface roughness of silicon on insulator (SOI) wafer.

#### **Scanning conditions**

System: NX-Wafer Scan Size: 1  $\mu$ m imes 1  $\mu$ m Scan Mode: Non-contact Scan Rate: 1 Hz

Cantilever: AC160TS (k=26 N/m, f=300 kHz) Pixel Size: 1024 × 256

#### **Scanning conditions**

System: NX20 Scan Size: Left 10 µm × 10 µm Right 10 µm × 10 µm Scan Mode: Non-contact Scan Rate: Left 2 Hz Right 2 Hz

#### True Non-contact<sup>™</sup> Mode



Height of position 2

Peak to valley: 1.4 nm RMS roughness: 0.15 nm

## Dendrimer



#### True Non-contact<sup>™</sup> Mode

In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.

## **Blood cell**



In this technique, the cantilever oscillates just above the surface as it scans. A precise, high-speed feedback loop prevents the cantilever tip from crashing into the surface, keeping the tip sharp and leaving the surface untouched. As the tip approaches the sample surface, the oscillation amplitude of the cantilever decreases. By using the feedback loop to correct for these amplitude deviations, one can generate an image of the surface topography.

#### Height (40 µm scan)



Peak to valley: 318.4 nm RMS roughness: 33.6 nm

## Height (2 µm scan)



Peak to valley: 209.6 nm RMS roughness: 22.6 nm

Peak to valley: 477.7 nm RMS roughness: 60.12 nm

#### Line profile



A late-stage erythroblast sample just above to expel its nucleus. The status of the erythroblast was confirmed by measuring the size of the nucleus, which has a width of 3 µm and peak to valley distance of 56.5 nm inside the nucleus. The Rq roughness inside the nucleus is 9 nm, which is flatter compared to the rest of the cell with a Rq of 25.8 nm.

#### **Scanning conditions**

System: NX10 Scan Size: Left 10 µm × 10 µm Right 2 µm × 2 µm

Scan Mode: Non-contact Scan Rate: Left 0.5 Hz Right 0.5 Hz

Cantilever: AC160TS (k=26 N/m, f=300 kHz) Pixel Size: Left 1024 × 1024 Right 1024 × 256

#### **Scanning conditions**

System: NX10 Scan Size: Left 40 µm × 40 µm Right 10 µm × 10 µm Scan Mode: Non-contact Scan Rate: Left 0.4 Hz Right 0.21 Hz

#### True Non-contact<sup>™</sup> Mode



## MLG-hBN, tBG+bBN (100 & 200 nm scan)



#### **CR-PFM**

PFM utilizes a lock-in amplifier to study the electrical properties and topography of a piezo sample surface in one single scan. Here, the AC voltage biased cantilever will introduce sample surface oscillation with same frequency. PFM signal can be enhanced by using the contact resonance frequency generated by the contact between the sample and the cantilever.

## Lithium niobate (LiNbO<sub>3</sub>, LN)



#### **PFM Amplitude on MLG-hBN**



PFM Amplitude on tBG-bBN





2D Moiré superlattice were measured using contact resonance PFM (CR-PFM).

MLG-hBN (Monolayer graphene on hexagonal boron nitride): 15 nm pitch honeycomb Moiré superlattice visible on the surface (Left)

tBG+bBN (Twisted bilayer graphene on boron nitride): 12 nm pitch honeycomb Moiré superlattice visible on the surface (Right)

• Sample courtesy:

- Qiong Ma, Boston College, US for tBG-bBN / David Goldhaber Gordon, Stanford University, US for MLG-hBN

#### **Scanning conditions**

System: NX10 Scan Size: Top 100 nm × 100 nm Bottom 200 nm × 200 nm Scan Mode: CR-PFM Scan Rate: Top left 1.5 Hz, right 3 Hz Bottom left 1 Hz, right 2 Hz

Cantilever: PPP-EFM (k=2.8 N/m, f=75 kHz) Pixel Size: Top left 512 × 256, right 256 × 256 Bottom left 256 × 256, right 256 × 512

#### **Scanning conditions**

System: FX40 Scan Size: 40 µm×40 µm Scan Mode: PFM Scan Rate: 0.4 Hz

#### **Piezoelectric Force Microscopy**

PFM utilizes a lock-in amplifier to study the electrical properties and topography of a piezo sample surface in one single scan. Here, the AC voltage biased cantilever will introduce sample surface oscillation with same frequency. The oscillation component of the PSPD signal is extracted by the lock-in amplifier, resulting in the PFM signal.

EL S



#### PFM Phase



# All-solid-state Li ion battery (ASS-LIB)

# All-solid-state Li ion battery (ASS-LIB)



## Kelvin Probe Force Microscopy

In Kelvin Probe Force Microscopy (KPFM), the AFM operates in non-contact mode while a conductive cantilever, oscillated at its fundamental resonant frequency, laterally scans over the sample surface. The resulting electrostatic signal provides information related to surface potential and the capacitance gradient. The topographic data is taken by controlling the force between the tip and the sample.



Schematic illustration of cross-sectional KPFM setup equipped with a source measure unit for (Cyclic voltammogram) CV operation of the ASS-LIB.



Schematic structure of the ASS-LIB. The region measured by KPFM is enclosed by a dashed square. Height



#### Surface potential during cyclic voltammetry operation





• Image courtesy: Dr. Nobuyuki Ishida, NIMS, Japan

#### **Scanning conditions**

System: XE-100 Scan Size: 10 µm × 4 µm Scan Mode: KPFM Scan Rate: 0.25 Hz



CV of the ASS-LIB during the KPFM measurement. The scan rate was 0.1 mV/s, and the voltage range was 0.2~2.7 V. COMMUNICATIONS CHEMISTRY | (2019) 2:140 | https://doi.org/10.1038/s42004-019-0245-x Dynamically visualizing battery reactions by operando Kelvin probe force microscopy Hideki Masuda, Kyosuke Matsushita, Daigo Ito, Daisuke Fujita & Nobuyuki Ishida

KPFM measurement on cathode composite region of ASS-LIB during CV operation. Surface potential images below are measured while the cell voltages were swept from 2.37 to 2.68 V.



• COMMUNICATIONS CHEMISTRY | (2019) 2:140

Cantilever: Multi75E-G (k=3 N/m, f=75 kHz) Pixel Size:  $1024 \times 64$ 

# PZT thin film



Sideband KPFM + CR-PFM after lithography

## Aluminium TX630 alloy









Height nm 200 100 -100 20 μm Peak to valley: 487.5 nm RMS roughness: 12.95 nm

#### **Scanning conditions**

System: FX40

Scan Size: 60 µm × 60 µm

Scan Mode: Sideband KPFM & CR-PFM after lithography - Lithography mode: Bias mode White/Black area +10 V/ -10 V bias

Scan Rate: Top 0.3 Hz Bottom 0.2 Hz Cantilever: PPP-NCSTAu (k=7.4 N/m, f=160kHz)

Pixel Size: All 1024 × 1024

**Scanning conditions** 

System: NX20 Scan Size: 50 µm × 50 µm Scan Mode: Sideband KPFM Scan Rate: 0.2 Hz

In Kelvin Probe Force Microscopy (KPFM), the AFM operates in non-contact mode while a conductive cantilever, oscillated at its fundamental resonant frequency, laterally scans over the sample surface. Sideband KPFM measures the surface potential of the sample using the signal that appears in the sideband of the resonant frequency. By using the force gradient, it measures the local interaction at the tip apex, not the average value acting on the entire lever.



Mean: 685 mV

Surface potential on Aluminium TX630 alloy prepared by semi-solid process.

## Cyanuric acid and melamine (CAM) on HOPG



Sideband KPFM + Tapping

# $F_{14}H_{20}$



Height



**Surface Potential** 





**Work Function** 

Phase



A melamine cyanuric network that forms a single layer on HOPG.

# Height

E 00 250 500 nm 750

Peak to valley: 5.45 nm RMS roughness: 1.12 nm

\* CAM: cyanuric acid (CA) and melamine (M)

#### **Scanning conditions**

System: NX20

Scan Size: 5 µm × 5 µm

Scan Mode: Top Tapping Bottom Sideband KPFM Scan Rate: All 0.25 Hz

Cantilever: Multi75E-G (k=3 N/m, f=75 kHz)

Pixel Size: All 2048 × 2048

#### **Scanning conditions**

System: FX40 Scan Size: 1 µm × 1 µm Scan Mode: Sideband KPFM Scan Rate: 0.5 Hz

In Kelvin Probe Force Microscopy (KPFM), the AFM operates in non-contact mode while a conductive cantilever, oscillated at its fundamental resonant frequency, laterally scans over the sample surface. Sideband KPFM measures the surface potential of the sample using the signal that appears in the sideband of the resonant frequency. By using the force gradient, it measures the local interaction at the tip apex, not the average value acting on the entire lever.



## **Surface Potential**

Peak to valley: 998.2 mV Mean: -108.6 mV

38

## **F**<sub>14</sub>**H**<sub>20</sub> (300 nm scan)



## Sideband KPFM

In Kelvin Probe Force Microscopy (KPFM), the AFM operates in non-contact mode while a conductive cantilever, oscillated at its fundamental resonant frequency, laterally scans over the sample surface. Sideband KPFM measures the surface potential of the sample using the signal that appears in the sideband of the resonant frequency. By using the force gradient, it measures the local interaction at the tip apex, not the average value acting on the entire lever.

# **Boron nitride on monolayer graphene**



#### Height



Peak to valley: 8.38 nm RMS roughness: 1.39 nm

#### **Surface Potential**



Peak to valley: 978.4 mV Mean: -155.5 mV

Height



Peak to valley: 13.8 nm RMS roughness: 0.89 nm

• Sample courtesy: Qiong Ma, Boston College, US

#### **Scanning conditions**

System: NX10	Scan Mode: Sideba
Scan Size: 1 $\mu m \times$ 1 $\mu m$	Scan Rate: 1 Hz

#### **Scanning conditions**

System: FX40 Scan Size: 300 nm × 300 nm Scan Mode: Sideband KPFM Scan Rate: 0.5 Hz

Cantilever: NSC36 Cr-Au C (k=0.6 N/m, f=65 kHz) Pixel Size: 512 × 256





## **Surface Potential**

Peak to valley: 613.9 mV Mean: -591.9 mV

Isolated square blocks were cut using voltage-biased contact mode (6 V<sub>Ac</sub>@17 kHz) and then imaged in sideband-KPFM.

## SRAM



#### **Conductive AFM**

The conductivity of the sample can be measured by performing a contact AFM scan with a conducting, biased tip. Regions of high conductivity on the sample surface alalow current to pass through easily, while regions of low conductivity will have a higher resistance. CP-AFM yields both the topography and the electrical properties of a sample surface.

1.2

Peak to valley: 2.2 nA Mean: -0.124 nA





#### Scanning conditions

System: NX20 Scan Size: 1.5 µm×1.5 µm Scan Mode: C-AFM Scan Rate: 1 Hz

Current measurement on SRAM with -1.5 V sample bias. P and N type of contact dot are well distinguished by IV spectroscopy measurements.

## Ta/ NiFe /Ta microman



#### Magnetic Force Microscopy

As much as EFM couples a topography scan with a simultaneous scan for electrical properties, Magnetic Force Microscopy (MFM) combines a topography scan with a concurrent scan for magnetic properties. MFM features a non-contact AFM scan to obtain the topography, and a scan farther from the surface to probe long-range magnetic force. In this magnetic force domain, deflections of the magnetized cantilever correspond.

## **PS-PMMA block copolymer**





Peak to valley: 107.6 nm RMS roughness: 14.6 nm

**MFM Phase** 









• Sample courtesy: Prof. Lew Wen Siang, NTU SPMS, Singapore

#### **Scanning conditions**

System: FX40 Scan Size: 10 µm × 10 µm Scan Mode: MFM Scan Rate: 0.5 Hz Cantilever: PPP-MFMR (k=2.8 N/m, f=75 kHz) Pixel Size: 512 × 256 Lift height: 50 nm

#### **Scanning conditions**

System: FX40 Scan Mode: Tapping, PinPoint nanomechanical mode

Scan Size: 6 µm × 6 µm

Scan Rate: Top Left/Middle 0.5 Hz Top Right/Bottom 0.05 Hz

Tapping + PinPoint<sup>™</sup> Nanomechanical Mode

Cantilever: Top Left/Middle PPP-NCSTAu (k=7.4 N/m, f=160 kHz) Top Right/Bottom BL-AC40 (k=0.09 N/m, f=110 kHz)

Pixel Size: Top Left/Middle 1024 × 512 Top Right/Bottom 512 × 256 43

# **Park SmartScan**<sup>™</sup>

Bringing the power and versatility of AFM technology to everyone





Park SmartScan<sup>™</sup> is a revolutionary operating software for Park AFMs that lets even inexperienced, untrained users produce high quality nanoscale imaging through three **simple clicks of a mouse** in auto mode, which rivals that made by experts using conventional techniques. SmartScan manual mode also provides all of the functions and tools necessary for more seasoned users to feel at home. This combination of extreme versatility, ease-of-use, and quality makes SmartScan the best AFM operating software available.

# **Park SmartLitho**<sup>™</sup>

The next generation nanolithography and nanomanipulation software combining powerful tools with an easy user interface



# **Park Systems**

## Dedicated to producing the most accurate and easiest to use AFMs

## General AFMs

Park Systems provides a range of popular AFMs for general research and industrial applications. Designed to be extremely versatile while still providing the accuracy and functionality necessary to do high quality work, our line of general AFMs offer researchers and engineers alike the ability to get extremely accurate results quickly and easily.



Park FX40 A Groundbreaking New Class of Atomic Force Microscope for Nanoscientific Research: The AutonomousAFM



Park NX10 The premiere chioce for nanotechnology research

#### • Materials Science • Failure Analysis

**Applications:** 

- Semiconductor Analysis
- Hard Disk Media Analysis



Park XE7 The most affordable research grade AFM with flexible sample handling



The most advanced high vacuum AFM for failure analysis and sensitive materials research

## **Industrial AFMs**

Park Systems is dedicated not just to advancing research, but industry as well. That's why our designers have worked to build a line of the most effective AFMs for FA engineers and industrial applications.



## Park NX-Hybrid WLI

The AFM and WLI technologies built into one seamless system



Park XE15 Power and versatility, brilliantly combined



Park NX20 The premiere choice for failure analysis



Park NX20 300 mm The leading automated nanometrology tool for 300 mm wafer measurement and analysis



Park NX-3DM

Automated industrial AFM for high-resolution 3D metrology



## Park NX-PTR

Fully automated AFM for accurate inline metrology of hard disk head sliders



## Park NX12 The most versatile AFM for analytical chemistry

## **Applications:**

- Failure Analysis
- Semiconductor Analysis
- Hard Disk Media Analysis



## Park NX-TSH

The automated Atomic Force Microscopy (AFM) system for ultra large and heavy flat panel displays at nanoscale



## Park NX-Wafer

Low noise, high throughput atomic force profiler with automatic defect review



## Park NX-HDM

Simply the best AFM for media & substrate manufacturing